

Test Update

IGH Somatic Hypermutation in B-CLL

UVM Medical Center

Effective 6/4/2026, IGH Somatic Hypermutation in B-Cell (Mayo BCLL) will be updated to clarify specimen requirements, and several result components will be added to allow for discrete results reporting. There will be no changes in ordering or billing.

Specimen Requirements

Current Specimen Requirements
Submit only 1 of the following:
Specimen Type: Whole blood
Container/Tube:
preferred: Lav Top (EDTA)
acceptable: Yellow Top (ACD)
Specimen Volume: 4 mL
Collection Instructions:
1. Invert several times to mix blood.
2. Send whole blood in original tube.
Do not aliquot.
3. Label specimen as blood.
Specimen Stability: Refrigerated/Ambient
Specimen Type: Bone marrow
Container/Tube:
preferred: Lav Top (EDTA)
acceptable: Yellow Top (ACD)
Collection Instructions:
1. Invert several times to mix bone marrow.
2. Send bone marrow in original tube.
Do not aliquot.
3. Label specimen as bone marrow.
Specimen Stability: Refrigerated/Ambient

New Specimen Requirements
Submit only 1 of the following:
Specimen Type: Whole blood
Container/Tube:
preferred: Lav Top (EDTA)
acceptable: Yellow Top (ACD)
Specimen Volume: 4 mL
Collection Instructions:
1. Invert several times to mix blood.
2. Send whole blood in original tube.
Do not aliquot.
3. Label specimen as blood.
Specimen Stability: Refrigerated (preferred) 7 days/Ambient 7 days
Additional Info: To ensure minimum volume and concentration of DNA are met, the requested volume must be submitted.
Testing may be canceled if DNA requirements are inadequate
Specimen Type: Bone marrow
Container/Tube:
preferred: Lav Top (EDTA)
acceptable: Yellow Top (ACD)
Specimen Volume: 2 mL
Collection Instructions:
1. Invert several times to mix bone marrow.
2. Send bone marrow in original tube.
Do not aliquot.
3. Label specimen as bone marrow.
Specimen Stability: Ambient (preferred) 7 days/ Refrigerated 7 days
Additional Info: To ensure minimum volume and concentration of DNA are met, the requested volume must be submitted.
Testing may be canceled if DNA requirements are inadequate.

UVM Medical Center

Pathology & Laboratory Medicine

Location

111 Colchester Ave.
Burlington, VT 05401

Lab Customer Service

802-847-5121
800-991-2799

Test Update

IGH Somatic Hypermutation B-CLL

UVM Medical Center

Test Build Information:

Orderable to be Updated	Epic Code	Atlas Code	Mayo Access ID	Order LOINC	CPT Code(s)
IGH Somatic Hypermutation in B-CLL	LAB16858	LAB16858	BCLL	50627-9	81263
Result Component(s)	Epic Code	Atlas Code	Mayo Access ID	Order LOINC	Notes
<i>Final Interpretation</i>	12301021539	M622044	622044	N/A	NEW
BCLL Result	12301015849	M39465	39465	N/A	
Specimen Type	12301015850	MMP005	MP005	31208-2	AOE
<i>Primer Set</i>	12301021540	M622055	622055	104271-2	NEW
<i>Primary V-gene</i>	12301021541	M622045	622045	104269-6	NEW
<i>Primary J-gene</i>	12301021542	M623097	623097	104270-4	NEW
<i>Primary Mutation Status %</i>	12301021543	M622046	622046	107552-2	NEW
<i>Primary Mutation Status</i>	12301021544	M622047	622047	107551-4	NEW
<i>Primary Rearrangement Status</i>	12301021545	M622048	622048	21747-1	NEW
<i>Primary Subset</i>	12301021546	M622053	622053	107550-6	NEW
<i>Secondary V-gene</i>	12301021547	M622049	622049	104269-6	NEW
<i>Secondary J-gene</i>	12301021548	M623098	623098	104270-4	NEW
<i>Secondary Mutation Status %</i>	12301021549	M622050	622050	107552-2	NEW
<i>Secondary Mutation Status</i>	12301021550	M622051	622051	107551-4	NEW
<i>Secondary Rearrangement Status</i>	12301021551	M622052	622052	21747-1	NEW
<i>Secondary Subset</i>	12301021552	M622054	622054	107550-6	NEW
Final Diagnosis	12301015851	M19674	19674	50398-7	

For questions regarding this change, contact UVMHC Lab Outreach at laboutreach@uvmhealth.org.