

Newly Defined Mayo Tests for UVMHC

Effective 12/1/22, the following newly defined Mayo tests will be available to order at UVMHC:

Lactoferrin, Quantitative, Stool

New Orderable	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code
Lactoferrin, QN, Stool	LAB17573	LAB17573	FLACQ	42924-1	83631
Result Component Name	Epic Code	Atlas Code	Mayo Test ID	Result LOINC	Notes
Lactoferrin, QN, Stool	12301019907	MFLACQ	FLACQ	42924-1	
Container	Specimen	Temperature	Collect Vol	Submit Vol	Stability
Sterile Stool Container	Fecal	Refrigerate	1 gm	1 gm	14 days
Collection Instructions					
Submit 1 gram fresh, unpreserved stool.					

PLA2R, Monitoring, ELISA, Serum

New Orderable	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code
PLA2R, Monitoring, ELISA, S.	LAB17558	LAB17558	PLA2M	73737-9	83520
Result Component Name	Epic Code	Atlas Code	Mayo Test ID	Result LOINC	Notes
PLA2R, Monitoring, ELISA, S.	12301019840	MPLA2M	PLA2M	73737-9	
Container	Specimen	Temperature	Collect Vol	Submit Vol	Stability
SST (preferred)	Serum	Refrigerate	2 mL	1 mL	14 days
Collection Information					
Red top tube is acceptable.					

Effective 12/19/22, the following newly defined Mayo test will be available to order at UVMHC:

Alzheimer's Disease Evaluation, CSF replaces LAB14512, Athena test #177:

New Orderable	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code
Alzheimer's Disease Evaluation, CSF	LAB17572	LAB17572	ADEVL	in process	83520x3
Result Component Name	Epic Code	Atlas Code	Mayo Test ID	Result LOINC	Notes
p-Tau/Abeta42	12301019881	MPTABR	PTABR	41027-4	
AD Interpretation	12301019882	MADINT	ADINT	69048-7	
Abeta42	12301019883	MAB42P	AB42P	33203-1	
Total-Tau	12301019884	MTTAUP	TTAUP	30160-6	
Phospho-Tau(181P)	12301019885	MPTAUP	PTAUP	72260-3	
Patient Preparation					
For 12 hours before specimen collection do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.					
Container	Specimen	Temperature	Collect Vol	Submit Vol	Stability
Blue Top SARSTEDT	CSF	Refrigerate	1.5 - 2.0 mL	1.5 - 2.0 mL	14 days
Collection Instructions					
1. Perform lumbar puncture and discard the first 1 to 2 mL of CSF.					
2. Collect CSF directly into a collection tube until the tube is at least 50% full.					
3. Send CSF specimen in original collection tube. Do not aliquot.					
Note: Polystyrene collection tubes are not acceptable. Exposure of CSF to polystyrene tubes may result in falsely low Abeta42 concentrations.					

PATHOLOGY & LABORATORY MEDICINE

111 Colchester Avenue | Mail Stop: 233MP1 | Burlington, Vermont 05401

PHONE LABORATORY CUSTOMER SERVICE

(802) 847-5121 | (800) 991-2799

Newly Defined Mayo Tests for UVMHC

Current Orderable to be Inactivated	Epic Code	Atlas Code	Athena Test ID	Order LOINC	CPT Code
ADMARK Phospho Tau/Total Tau/A Beta42, Analysis and Interp, CSF	LAB14512	N/A	#177	N/A	83520x3

For questions or concerns, please contact UVMHC Lab Outreach at laboutreach@uvmhealth.org.

PATHOLOGY & LABORATORY MEDICINE

111 Colchester Avenue | Mail Stop: 233MP1 | Burlington, Vermont 05401

PHONE LABORATORY CUSTOMER SERVICE

(802) 847-5121 | (800) 991-2799