

TEST UPDATE

Creutzfeldt-Jakob Disease Evaluation, CSF

UVMMMC

Effective 12/11/24, the Mayo test Creutzfeldt-Jakob Disease Evaluation, CSF (Mayo CJDE) will be available to order using a distinct orderable code. This has been defined for ease of ordering and to assist practitioners in specimen collection.

Test Build Information:

New Orderable	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code(s)
Creutzfeldt-Jakob Disease Evaluation, CSF	LAB17725	LAB17725	CJDE	97502-9	84999, 83520x2
Result Component(s):	Epic Code	Atlas Code	Mayo Test ID	Result LOINC	Notes
CJD Eval Interp, CSF	12301020441	M620375	620375	69048-7	
RT-QuIC Prion, CSF	12301020443	M620307	620307	101662-5	
t-Tau/p-Tau	12301020444	MCTTPT	CTTPT	101752-4	
Total-Tau	12301020445	MCTTAU	CTTAU	30160-6	
Phospho-Tau(181P)	12301020446	MCPTAU	CPTAU	72260-3	
Specimen Requirements:					
Container	Specimen	Temperature	Collect Vol	Submit Vol	Stability
2 Sardstedt CSF False Bottom Tubes*	CSF	Frozen	2.5 mL ea.	1.5 mL ea.	28 days
Collection Information:					
*Included in CJD/RPD Evaluation Kit (T966)					
1. Perform lumbar puncture and discard the first 1 to 2 mL of CSF.					
2. Collect 2 tubes of CSF directly into an acceptable collection tube until the tube is at least 50% full.					
3. Send CSF specimen in original collection tube. Do not aliquot.					

For questions regarding this change, contact UVMMMC Lab Outreach at laboutreach@uvmhealth.org.

PATHOLOGY & LABORATORY MEDICINE

111 Colchester Avenue | Mail Stop: 233MP1 | Burlington, Vermont 05401

PHONE LABORATORY CUSTOMER SERVICE

(802) 847-5121 | (800) 991-2799