

TEST UPDATE

IA-2 Autoantibodies Testing Change

UVMMC

Effective 8/18/25, IA-2 Autoantibodies (Mayo FIAAB) will replace Islet Antigen 2 (IA-2) Antibody, Serum (Mayo IA2). Islet Antigen 2 (IA-2) Antibody, Serum is being discontinued as part of ongoing efforts to enhance testing efficiencies and reduce the use of radioimmunoassays.

Test Build Information:

New Orderable	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code(s)
IA-2 Autoantibodies	LAB17823	LAB17823	FIAAB	31209-0	86341
Result Component(s):	Epic Code	Atlas Code	Mayo Test ID	Result LOINC	Notes
IA-2 Autoantibodies	12301020961	MFIAAB	FIAAB	31209-0	
Specimen Requirements:					
Container	Specimen	Temperature	Collect Vol	Submit Vol	Stability
SST	Serum	Frozen	2 mL	1 mL	107 days
Collection Information:					
Red Top tube is acceptable.					
Minimum volume: 0.5 mL (doesn't allow for repeat testing)					
Freeze immediately after separation and send frozen.					
Refrigerated and Ambient temps are acceptable, but specimens are only stable for 72 hours.					

Orderable to be Discontinued	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code(s)
Islet Antigen 2 (IA-2) Antibody, S.	LAB12112	IA2	IA2	81155-4	86341

For questions regarding these changes, contact UVMMC Lab Outreach at laboutreach@uvmhealth.org.

PATHOLOGY & LABORATORY MEDICINE

111 Colchester Avenue | Mail Stop: 233MP1 | Burlington, Vermont 05401

PHONE LABORATORY CUSTOMER SERVICE

(802) 847-5121 | (800) 991-2799

THE
University of Vermont
MEDICAL CENTER