

Test Update

Lactoferrin, Feces Test Changes

UVM Health

Effective 1/15/26, the new Mayo test, Lactoferrin, Feces (Mayo LCTF) will replace Lactoferrin, Fecal by ELISA (Mayo FLACF) and Lactoferrin, Quantitative, Stool (Mayo FLACQ), both of which will be obsolete.

Test Build Information:

New Orderable	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code(s)
Lactoferrin, Feces	LAB17895	LAB17895	LCTF	in process	83631
Result Component(s):	Epic Code	Atlas Code	Mayo Test ID	Result LOINC	Notes
Lactoferrin, Feces	12301021414	MLCTF	LCTF	in process	
Specimen Requirements:					
Container	Specimen	Temperature	Collect/Submit Vol	Min Vol	Stability
Stool Container	Feces	Frozen	5 grams	1 gram	7 days
Collection Instructions:					
Ambient and Refrigerated specimens are acceptable. Stability for these is 72 hours.					
Refrigerated or thawed specimens received more than 72 hours after collection will be rejected.					
Testing cannot be added-on to a previously collected specimen.					

Orderable to be Discontinued	Epic Code	Atlas Code	Mayo Test ID	Order LOINC	CPT Code(s)
Lactoferrin, Fecal by ELISA	LAB15576	LAB15576	FLACF	40703-1	83630
Lactoferrin, Quantitative, Stool	LAB17573	LAB17573	FLACQ	42924-1	83631

For questions regarding these changes, please contact UVMHC Lab Outreach at laboutreach@uvmhealth.org.